# **Ordering Information**

Part Number

Description

#### Intelect Mobile Laser

Intelect Mobile Laser

#### Standard Accessories

27525 Laser Protection Glasses (2x)

28063 User Manual on CD

For Laser Applicators, please see page 18

# Optional Accessories

27904K	Remote Interlock Kit
27467	Intelect Mobile Carrying Bag
27478	Battery Pack
2775	Therapy System Cart
2883	Cart Adapter

# **Technical Specifications**

Mains Power: 120-240 VAC, 50 60 Hz

Weight: 2.3 kg

Dimensions\*: 33 X 29 X 16.3 CM

Electrical

Class 1, Type BF Safety Class:

Safety Tests: IEC 60601-1, IEC 60825-1, IEC 60601-2-22

# Intelect® Mobile Laser





#### Features

- Class 3B Low Level Laser Therapy
- Independent control over all parameters
- Continuous and pulsed operation available
- Fully variable pulsed frequencies
- Real time feedback of energy delivered
- Display of energy in choice of joules or joules/ cm²
- Helpful menu of Clinical Indications specific to the laser applicator
- 15 user-defined memory slots on board

- Built in acupuncture point finder
- 17 different applicators available, including LED, laser diodes, and clusters
- Innovative design allows for cart, table top, wall or mobile use
- Backlit LCD display for better visibility in low light settings
- Improved ▲  $\blacktriangledown$  button for easier navigation
- 3 year warranty (probes and accessories 1 year)



<sup>\*</sup> Dimensions are listed as L x W x H